

NIC Nano hydro

REFINE THE PATH

IN TIGHTEST LESIONS

 **0.85 mm**

LOW-PROFILE
NON-COMPLIANT PTCA
DILATATION CATHETER

High RBP 21 atm

Lubricious hydrophilic coating

Marker

Single Pt/Ir proximal marker

Tip


Unique tip design to facilitate navigation

**MADE FOR TIGHTEST
LESIONS AND CTO**

NIC Nano hydro

0.85 MM LOW-PROFILE NON-COMPLIANT PTCA DILATATION CATHETER

Balloon Diameter, mm	Balloon Length, mm	Art. Number	UDI-DI-Number (GTIN)
0.85	6	085-006-134	7640142573029
0.85	10	085-010-134	7640142573036
0.85	15	085-015-134	7640142573043



TECHNICAL SPECIFICATIONS

Nominal Pressure	10 atm
Rated Burst Pressure (RBP)*	21 atm
Guiding Catheter Compatibility	5F
Guide Wire Compatibility	0.014" (0.36 mm)
Marker position	Single Pt/Ir proximal marker
Catheter Usable Length	140 cm
Catheter Coating	Hydrophilic

*Do not exceed RBP

NIC Nano hydro is not approved for sale in the USA.

SIS MEDICAL

Hungerbühlstrasse 12a
CH-8500 Frauenfeld
Switzerland
+41 52 245 09 90
sales@sis-medical.com
www.sis-medical.com

Illustrations are not to be recognized as engineering drawings.

GERMANY
Distributed by: SIS Medical GmbH
Kirchstrasse 36, 72127, Kusterdingen, Germany
+49 7127-799 8830
orders-de@sis-medical.com