

# EasyT

FOR EXCELLENT LESION PREPARATION, EVEN IN MOST COMPLEX LESIONS

## WIDE VARIETY OF SIZES

### Low crossing profile

Crossing profile 0.035 inch ( $\varnothing$  3.0 mm)

### Markers

Dual Pt/Ir markers for all balloon sizes

### Tip

Unique tip design  
to facilitate navigation

### Folding technology

Tight trifold\* to enhance crossability

MADE WITH **SWISS PRECISION**



**SWISS MADE**

[www.sis-medical.com](http://www.sis-medical.com)

\*EasyT 1.5 mm diameter is bi-fold

# EasyT

## HIGH-PRESSURE PTCA DILATATION CATHETER

### BALLOON DIAMETERS

Art. Numbers	1.50 mm	2.00 mm	2.50 mm	3.00 mm	3.50 mm	4.00 mm	4.50 mm
<b>6 mm</b>	150-006-134	200-006-134	250-006-134	300-006-134	350-006-134	400-006-134	450-006-134
<b>10 mm</b>	150-010-134	200-010-134	250-010-134	300-010-134	350-010-134	400-010-134	450-010-134
<b>15 mm</b>	150-015-134	200-015-134	250-015-134	300-015-134	350-015-134	400-015-134	450-015-134
<b>20 mm</b>	150-020-134	200-020-134	250-020-134	300-020-134	350-020-134	400-020-134	450-020-134
<b>25 mm</b>		200-025-134	250-025-134	300-025-134	350-025-134	400-025-134	450-025-134
<b>30 mm</b>		200-030-134	250-030-134	300-030-134	350-030-134	400-030-134	450-030-134

## BALLOON COMPLIANCE CHART

### BALLOON DIAMETERS

		1.50 mm	2.00 mm	2.50 mm	3.00 mm	3.50 mm	4.00 mm	4.50 mm		
PRESSURE atm	<b>10</b> NP	1.55	2.02	2.45						
	<b>12</b>	1.59	2.11	2.54	3.01	3.43	3.95	4.45	NP	
	<b>14</b>	1.63	2.17	2.61	3.12	3.55	4.11	4.61		
	<b>16</b>	1.66	2.21	2.66	3.20	3.65	4.23	4.73		
	<b>18</b>	1.68	2.24	2.71	3.26	3.72	4.32	4.82		
	<b>20</b>	1.70	2.27	2.76	3.31	3.79	4.39	4.89		
	<b>21</b> RBP	1.73	2.29	2.78	3.35	3.83	4.44	4.93	RBP	

## TECHNICAL SPECIFICATIONS

	1.5 - 2.5 mm	3.0 - 4.5 mm
Nominal Pressure (NP)	10 atm	12 atm
Rated Burst Pressure (RBP)**	21 atm	21 atm
Guiding Catheter Compatibility	5F	5F
Guide Wire Compatibility	0.014" (0.36 mm)	0.014" (0.36 mm)
Marker	Two Pt/Ir markers	Dual Pt/Ir markers
Catheter Usable Length	140 cm	140 cm
Catheter Coating	Hydrophilic	Hydrophobic

SIS MEDICAL

Hungerbuelstrasse 12a  
CH-8500 Frauenfeld  
Switzerland  
+41 52 245 09 90  
sales@sis-medical.com  
www.sis-medical.com

GERMANY  
Distributed by: SIS Medical GmbH  
Kirchstrasse 36, 72127, Kusterdingen, Germany  
+49 7127-799 8830  
orders-de@sis-medical.com

Manufactured by SIS Medical AG.  
All rights reserved.

\*\*Do not exceed RBP

EasyT is not approved for sale in the USA.

Illustrations are not to be recognized as engineering drawings.